

GB

- Tight economy
- Scarce resources

Pressure Ulcers

Try Levabo's product range of inflatable one-patient cushions

- Hygienic
- Comfortable
- Easy
- Economic



LEVABO[®]
medico cushions



Cushions for pressure relief of heels and legs



Name	Size in cm.	Item no.	Quantity	Application	About the cushion
Heel up - Short	18 x 32 x 25	72001	12 pcs		Heel protector for mobile and seated users, typically wheelchair. Can be used for both prevention and treatment of pressure sores.
Heel up - Medium	18 x 39 x 25	72007	12 pcs		Heel protector useful for both prevention and treatment of pressure sores. Typically for users with short legs. Most useful for bed-bound patients. It supports the front foot and minimizes risk of foot drop.
Heel up - Long	18 x 48 x 25	72002	12 pcs		Heel protector useful for both prevention and mitigation of pressure sores. Most useful for bed-bound patients. It supports the front foot and minimizes the risk of foot drop.
Heel up - Max	18 x 50 x 30	72018	12 pcs		Heel protector useful for both prevention and mitigation of pressure sores. Most useful for bed-bound patients and for users with a long foot. It has an extra air chamber under the foot, which provides an optimal support of the front foot.



All up - Multi Small	50 x 25 x 8	72003	12 pcs		Positioning cushion with 2 chambers. Multifunction. Can be used for both supine and sitting positions. (between legs, back, very good support for the arm and between the knees).
All up - Multi Medium	50 x 27 x 8	72004	12 pcs		Positioning cushion with 3 chambers. Multifunction. Especially useful for a supine position (between legs, back), and possibly as support under the feet.
All up - Multi Large	50 x 50 x 8	72005	12 pcs		Positioning cushion 4 chambers. Multifunction. With the possibility of being folded/rolled up. (Under the leg and arm; comes with extra velcro strap).
All up - Multi Max	50 x 50 x 20	72012	6 pcs		Positioning cushion, particularly useful for hip fracture patients and when the lower limb need to be elevated above heart level.
All up - Wing	50 x 85 x 8	72006	6 pcs		Relief cushion for lower limbs, used for bed-bound patients. The cushion is easily fixed to the mattress.

Accessories

Name	Size in cm.	Item no.
Heel up - Inlay X 10	35 x 50	76008
All up - Cover Medium X 5	50 x 50	76012
All up - Cover Large X 5	50 x 100	76013





Awards


Levabo was awarded The Innovation Elder Care Award 2015 (Innovationspreis Altenpflege 2015) for the product series Heel up® on Germany's largest elder care conference.

The panel of judges especially emphasized how easy applicable the products are, as the inflator indicator Stop Up® ensures correct inflation without the user being required to have extensive training on how to inflate the cushions.

Additionally, the products were praised for their compact size not taking up space in the depots, and them being easy to transport to patients, particularly and advantage for the domiciliary care.



Remember the exposed areas



The first sign of a pressure sores can develop within a very short time, 2 – 4 hours !

Pressure sores typically develop on the exposed areas of the body, where the bones are only protected by a thin layer of skin.

With Levabo's flexible one-patient product range it is possible to make a quick and effective positioning of the patient before the exposed areas are affected.

Cushions for pressure relief in the seated position



Name	Size in cm.	Item no.	Quantity	Application	About the cushion
All up - Seat	42 x 42 x 5	72008	12 pcs		Seat cushion, extra strong. Pressure distribution for the mitigation of pressure sores. Standard wheelchair size. Max 150 kg.
All up - Seat Large	52 x 52 x 5	72019	12 pcs		Seat cushion, extra strong. Pressure distribution for the prevention and mitigation of pressure sores. Extra length and width. Max 150 kg.
All up - Seat Max (bariatrisk)	52 x 62 x 5	72020	12 pcs		Seat cushion, extra strong. Bariatric cushion for user weights up to 250 kg. Pressure distribution for the prevention and mitigation of pressure sores. Extra length and width.
All up - Seat Twin	48 x 52 x 5	72021	12 pcs		Seat cushion, extra strong. Suitable for hemiplegic patients with unstable posture. The cushion is split in two to adjust for individual air pressure. For the prevention and relief of pressure sores. Max 150 kg.
All up - Seat Tailbone	42 x 42 x 5	72022	12 pcs		Special sitting cushion, extra strong. To apply for full relief of the rectum, typically for radiation damage, a sensitive tailbone or great discomfort caused by hemorrhoids. Max 150 kg.
All up - Back & Seat	27 x 45 x 5	72009	12 pcs		Sitting cushion and back cushion. With advantage, usefull for both areas. Good back support, easy cushion for correct seating. Max 150 kg.
All up - Donut	45 x 45 x 8	72017	12 pcs		Seat cushion, extra strong. Annular (ring shaped) with no hole in the middle, but with optimum pressure to the relief and prevention of pressure sores. Max 150 kg.
All up - Bath	50 x 50 x 8	72013	12 pcs		Special seat cushion, extra strong. For bath and toileting chairs to offer a better distribution of pressure for prolonged use. Can be fastened to the chair with tape. After bathing, the cushion can be dried and reused. Max 150 kg.

Prevention and treatment

Cushions for special positioning

All up - Belly	65 x 75 x 15	72023	6 pcs		Special positioning cushions with particular application to the radiology department to lift the chest at special examinations.
All up - Kidney	50 x 20 x 20	72014	12 pcs		Special positioning cushion for special for a Radiology Department for seating kidney/pelvic area.

Levabo Emergency Package

Our products really make a difference.

The innovative positioning products from Levabo[®] are developed in collaboration with leading wound care nurses and are produced in Denmark. The products have a series of advantages compared to conventional positioning products.


With the use of Levabo[®]'s positioning cushions you always have the needed cushion at hand for effective prevention of pressure sores. The products are one-patient, inflatable positioning cushions, which can be disposed after use. The products take up very little space, both useful for the depot, but also for easy transportation of Levabo[®]'s emergency package, developed for that purpose.

These advantages are especially useful for the domiciliary care. As a domiciliary care nurse you can initiate effective prevention of pressures sores at an early stage with Levabo[®]'s products, easier and cheaper than with conventional positioning products. Already at suspicions of the development of pressure sores, the domiciliary nurse can immediately position the patient optimally, and thereby avoid unnecessary discomfort for the patient, as well as save the health care system for potential treatment expenses.

The products are multifunctional and can with advantage be used in several different positions. The one-patient use ensures high hygiene.



Name	Size in cm.	Item no.
dcubi kit bag	40 x 31 x 15	76007








Prevention of pressure sores

Accessories



Name	Size in cm.	Item no.
All up - Seat Cover X 5	42 x 42 x 5	76009
All up - Seat Large Cover X 5	52 x 52 x 5	76010
All up - Sheet Medium X 5	50 x 50	76012
All up - Sheet Large X 5	50 x 100	76013

Cushions for pressure relief of body and arms

Name	Size in cm.	Item no.	Quantity	Application	About the cushion
All up - Multi Small 	50 x 25 x 8	72003	12 pcs		Positioning cushion with 2 chambers. Multifunction. Can be used for both supine and sitting positions. (between legs, back, very good support for the arm and between the knees).
All up - Multi Medium 	50 x 27 x 8	72004	12 pcs		Positioning cushion with 3 chambers. Multifunction. Especially useful for a supine position (between legs, back), and possibly as support under the feet.
All up - Multi Large 	50 x 50 x 8	72005	12 pcs		Positioning cushion 4 chambers. Multifunction. With the possibility of being folded/rolled up. (Under the leg and arm; comes with extra velcro strap).
All up - Multi Max 	50 x 50 x 20	72012	6 pcs		Positioning cushion, particularly useful for hip fracture patients and when the lower limb need to be elevated above heart level.
All up - Universal 	50 x 50 x 8	72010	12 pcs		Positioning cushion, typically for the back, sacrum and tailbone in adjustable bed and the lateral positions between the legs. Can be folded if a higher cushion is desired.
All up - Universal Long 	55 x 70 x 8	72015	12 pcs		Positioning cushion, typically for the lateral positions between the legs, longer model which protects the knees and malleolus. Can be folded in half if a higher cushion is desired.
All up - Universal Slim 	45 x 75 x 8	72016	12 pcs		Positioning cushion, typically for lateral positions between the legs, longer and narrower model which ensures a distance between the knees in the lateral position. Can be folded if a higher cushion is desired.
All up - Head 	28 x 30 x 5	72011	12 pcs		Cushion, particularly useful for the preventing and treating of pressure sores on the ears and neck, and is also useful for coma patients.
All up - Wing 	50 x 85 x 8	72006	6 pcs		Relief cushion for lower limbs, used for bed-bound patients. The cushion is easily fixed to the mattress.
All up - Wedge 30 deg. 	70 x 40 x 15	72024	6 pcs		Positioning cushion for the lateral position. Ensures a stable position at 30°, and the wedge shape is fixed by a flap that is pushed underneath the patient.
All up - Lap 	80 x 60 x 10	72025	6 pcs		Positioning cushion, crescent. Extra height. Advantageously used for wheelchair users for relief and support of arms.
All up - Lap Max 	80 x 60 x 20	72026	6 pcs		Positioning cushion, crescent. Advantageously used for wheelchair users for relief and support of arms.
All up - Arm Rest 	55 x 25 x 8	72027	12 pcs		Positioning cushion for relief of arms.
All up - Elbow 	30 x 20 x 10	72028	12 pcs		Positioning cushion for relief of the elbow.



Levabo's products are space saving in the depot.



Accessories

Name	Size in cm.	Item no.
All up - Strap X 1 (Velcro)	5.5 x 200	76011
All up - Cover Medium X 5	50 x 50	76012
All up - Cover Large X 5	50 x 100	76013

Levabo Pumps


All Levabo[®] products are delivered with an inflation adapter set, which is included in the packaging of the individual cushion.

Name	Size in cm.	Item no.
Levabo Manual Pump	29 x 6 x 6	76004
Levabo Battery Pump	26 x 5 x 5	76001
Levabo 240V Pump	24 x 13 x 11	76002
Levabo inflate adapter set X 10		76005



Always correct inflation

With Levabo®'s Stop Up® patent system you are always sure to make a correct inflation. The Stop Up® strip indicates when the right amount of air is in the cushion. Check out the video at: www.levabo.com



Levabo®'s Stop Up® system with easily to understandable pictograms.

- Single-patient use ensures increased hygiene.
- Heel up® is designed for the treatment and prevention of bedsores in the heel region.
- All up® was designed for the prevention and treatment of bedsores on the elbows, shoulders, back, etc.
- Easy & simple to use.
- Good comfort
- Soft, skin-friendly, heat dissipation and moisture-absorbent material.
- Space saving in the storage area.
- Delivered in clean, sealed and flat packaging, including an inflation adapter set.
- Tested and approved according to ISO 10993-10 skin irritation and hypersensitivity.
- Available in several sizes.

Levabo Guide

For your iPhone, iPad & Android



Get the free Levabo app on:

<http://www.levabo.dk/apps.html>

Or scan the QR code:



Dealer



www.levabo.dk
www.facebook.com/LevaboAps



LEVABO®
medico cushions

Levabo ApS
Højlundsparken 49
8355 Solbjerg
Tel.: +45 20 70 09 30
info@levabo.dk
www.levabo.com